

# Hemorrhoids

## Treatment

Hemorrhoids which are not treated properly can gradually develop up to degree 3 or 4, which could result into an urgent need of a painful operation. With the HPR45i probe, all types of hemorrhoids can be treated easily and effectively!


This new and improved method is easy to use and very effective on hemorrhoids degree 1-4. Treatments can be done in one or two sessions of 15-20 minutes under local anesthesia. Patients hardly experience any pain after treatment so minimal postoperative care is required and a swift return to normal daily activities is possible.



## 4 Step procedure



1. Place the proctoscope in the rectum of the patient to see the hemorrhoid.



2. Apply local anaesthesia (Lignocaine 1%, 3 to 6 ml) directly in the hemorrhoid to be treated.



3. Insert the HPR45i probe in the hemorrhoid and apply during 5 to 10 beeps 25-watt.



4. Repeat until the hemorrhoid is completely coagulated.

## Clinical trial

Effectiveness and pain of hemorrhoid treatments stage 3 and 4 with the EVRF<sup>®</sup> on 300 patients.

Kryocentrum Poland,  
Dr. Klos.

# EVRF<sup>®</sup>

ENDOVENOUS RADIO FREQUENCY

## Needles and catheters

### The advantages

#### For the patient

- Minimal discomfort during the treatment
- No or minimal bleeding after the treatment
- Quick return to daily activities
- Only local anesthesia
- Walk-in walk-out daycare treatment


#### For the doctor

- Only local anesthesia
- Easy to learn procedure and easy to apply
- Immediate clinical results
- No post-operative care
- Cost effective and efficient treatment
- Day clinic or private practice

### Technical

#### specifications HPR45i

- |            |                   |
|------------|-------------------|
| • Diameter | 1.9 mm            |
| • Length   | 200 mm            |
| • Weight   | 34 gr             |
| • Coating  | PTFE              |
| • Tip      | Sharp Metal 50 mm |



F Care Systems NV  
Kontichsesteenweg 54  
2630 Aartselaar • Belgium  
Tel: +32 3 451 51 45  
Fax: +32 3 451 51 39  
info@fcaresystems.com  
[www.fcaresystems.com](http://www.fcaresystems.com)

