

Varicose veins (from I to 5 mm)

Treatment

Reticular, collateral and perforating veins can be treated with F Care Systems' CR12i and CR30i catheters. The catheters are extremely flexible, so that they follow the direction of the vein easily. Smooth insertion is ensured with the advanced coating material around the catheter. The non-insulated tip transmits the high frequency signal to the vein wall. This causes the vein to coagulate and eventually disappear.

The catheter is manipulated with a sterilized handset. Scrolling the handset's wheel will move the catheter in or out of the vein.



CR30i catheter with handset

4 Step procedure







2. Replace the introduction needle with the catheter using the sterilized handset.



3. Move the catheter up the vein by scrolling the handset wheel.



 Retract the catheter while pushing the treadle to transmit an RF signal in the vein. The vein will coagulate and disappear.





Needles and catheters

The advantages

For the patient

- No pain during and after treatment
- No pigmentation or allergic skin reaction
- · No dermatological reactions due to chemical compound
- High esthetic outcome

For the doctor

- Easy to use
- Permanent occlusion
- All body parts
- Semi-sterile environment (private practice)
- · No chemical agent in the body
- · Collateral, peripheral and perforating veins can be treated

Technical specifications

CR12i

Diameter catheterLength27 mm

• Handset 140 mm x 55 mm x 15 mm

Weight handset 40 grCoating PTFETip Metal

CR30i

Diameter catheterLength27 mm

Handset 140 mm x 55 mm x 15 mm

Weight handset 40 grCoating PTFETip Metal

F Care Systems NV
Kontichsesteenweg 54
2630 Aartselaar • Belgium
Tel: +32 3 451 51 45
Fax: +32 3 451 51 39
info@fcaresystems.com
www.fcaresystems.com

