

6 ANNUAL CONTROL FORMS

The **EVRF** must have an annual control to be able to claim with the maintenance of its MARKING EC. This control makes it possible to check all the functions of the apparatus and, particularly, the safety measures.

This chapter relates to the notebook of maintenance of the EVRF.

Control of end of guarantee	
Date of control to carry out:	Date of control :
Checks carried out: <ul style="list-style-type: none"> <input type="checkbox"/> - Visual monitoring and cleaning of EVRF. <input type="checkbox"/> - Security's control. <input type="checkbox"/> - Control functions of orders. <input type="checkbox"/> - Measuring of the output signals. <input type="checkbox"/> - Control of the foot pedal. <input type="checkbox"/> - Control pen carries needle. 	
Changed components: <ul style="list-style-type: none"> <input type="checkbox"/> - Pen carried needle. <input type="checkbox"/> - Other components : 	
Name and visa of the controller	Name and visa of the customer: